

TECHNICAL PROTOCOL  
FOR  
**294-Cre/ 294-FLP**  
*E. coli* strain  
(version 2.0)

Gene Bridges GmbH  
Pfortenhauerstr. 108  
01307 Dresden, Germany  
Tel + 49 (0)351 210 2552  
Fax + 49 (0)351 210 1619  
Email: [products@genebridges.com](mailto:products@genebridges.com)  
[www.genebridges.com](http://www.genebridges.com)

## CONTENTS

### 1 Eppendorf tubes + manual

1. 294-Cre or 294-FLP: Glycerol stock of the *E. coli* strain (500µl, 25% glycerol)
2. This manual

### Store tube at -80°C

### Please read

The products listed in this manual are for research purposes only. They are not designed for diagnostic or therapeutic use in humans, animals or plants. The Red<sup>®</sup>/ET<sup>®</sup> recombination technology is the intellectual property of Gene Bridges GmbH.

### Conditions of use

3.1 Purchaser will not manufacture, copy, reproduce, transmit, distribute, sell, lease, transfer, or improve upon the MATERIALS without prior written consent from GENE BRIDGES.

3.2 All MATERIALS relating Technologies shall be purchased from GENE BRIDGES or its authorized distributors. Use of any of the stated products from a source other than GENE BRIDGES will exempt GENE BRIDGES from any and all liabilities and warranties.

3.3 All MATERIALS purchased by research organizations, universities and other non-profit organizations may not be used for any commercial purpose. These MATERIALS are to be used for research purposes only. The MATERIALS may not be used to provide a commercial or non-commercial service, of any kind.

3.4 A purchase of MATERIALS by a private consumer is neither intended nor permitted.

**Short Description:**

*E.coli* strains 294-Cre and 294-FLP express either Cre- or FLP-recombinase. Plasmids containing authentic recognition targets for either recombinase (loxP or FRTs) are recombined when propagated in the appropriate strain. 294-Cre and 294-FLP thus provide a simple test for the recombination competence of constructs that are designed for use in Cre- or FLP-mediated genomic manipulations.

In both strains the recombinases are under the control of the  $\lambda P_R$ -promoter which limits expression at 23°C so that no FLP, and little Cre recombination is observed in cultures grown at this temperature. Plasmids with recombinase recognition targets are completely recombined at 37°C in the appropriate Cre or FLP strain.

*E.coli* strains 294-Cre and 294-FLP are based on the MM294 wild-type strain.

Genotype MM294: endA, thiA, hsdR17, supE44

## Literature:

- Buchholz F., Angrand P.O. and Stewart A.F. 1996: A simple assay to determine the functionality of Cre or FLP recombination targets in genomic manipulation constructs. *Nucleic Acids Research* 24, 3118-3119.
- Morschhäuser J., Michel S. and Staib P. 1999: Sequential gene disruption in *Candida albicans* by FLP-mediated site-specific recombination. *Molecular Microbiology* 32, 547-556.
- Lowe L.A., Yamada S. and Kuehn M. R. 2001: Genetic dissection of nodal function in patterning the mouse embryo. *Development* 128, 1831-1843.
- Reizis B. and Leder P. 2001: The upstream enhancer is necessary and sufficient for the expression of Pre-T cell receptor  $\alpha$  gene in immature T lymphocytes. *J. Exp. Med.* 194, 979-990.
- Langer S. J., Ghafoori A. P., Byrd M. and Leinwand L. 2002: A genetic screen identifies novel non-compatible loxP sites. *Nucleic Acids Research* 30, 3067-3077.
- Schnütgen F., Doerflinger N., Calleja C., Wendling O., Chambon P. and Ghselinck N.B. 2003: A direct strategy for monitoring Cre-mediated recombination at the cellular level in the mouse. *Nature Biotechnology* 21, 562-565
- Seibler J., Zevnik B., Kuter-Luks B., Andreas S., Kern H., Hennek T., Rode A., Heimann C., Faust N., Kauselmann G., Schoor M., Jaenisch R., Rajewsky., Kuhn R. and Schwenk F. 2003: Rapid generation of inducible mouse mutants.
- Zhang X.-M. and Huang J.-D. 2003: Combination of overlapping bacterial artificial chromosomes by a two-step recombinogenic engineering method. *Nucleic Acids Research* 31, e81.

## Technical Protocol

1. Streak 1µl of the cells on an LB agar plate (without antibiotics).
2. Incubate at 30°C overnight.
3. Pick several colonies and grow them in LB medium (without antibiotics). Incubate them at 30°C.
4. Prepare competent cells by any conventional method (electro-competent cells are preferable).
5. Transform the plasmid carrying loxP or FRT sites into 294-Cre or 294-Flp competent cells.
6. After transformation, add 1 ml of LB medium.
7. Incubate at 30°C for 1.5 hours.
8. Plate on the appropriate selection plates and incubate at 30°C for more than 24 hrs.
9. Pick colonies from the plates and isolate plasmid DNA following any of the conventional methods.
10. Re-suspend the DNA pellet in 10µl of 10mM Tris-HCl, pH8.0.
11. Transform 1-2 µl of the plasmid DNA into another *E.coli* strain (DH10B, Top10, DH5α ...).
12. Isolate plasmid DNA and check the successful recombination event by digestion.

Since 294-Cre and 294-FLP are both based on the MM294 wild-type *E.coli* strain, the plasmid DNA prepared from them is not of high quality. Therefore, after recombination, the plasmid should be transformed into another strain (e.g. DH10B, Top10, DH5α...).