



Product: Mouse embryonic stem cells 17.2

Catalog #: MES1429 Lot#:

Cell storage

-140°C indefinitely

Product Description

- Quantity: 1.0ml (approximately 5×10^6 cells/cryovial)
- Cell passage number: 18
- Strain: BALB/c x 129/sv (M)

Quality Control

- All cell lines were determined to be mycoplasma negative.
- Post-cryostorage viability was determined to be approximately 80%.
- Germline transmissibility tested at passage number 8.

Usage recommendations

We recommend thawing 1 cryovial per 1 well of a 6 well plate (9.5cm²)