



---

**Product: Mouse NH-MEF Large**

Catalog #: MES3950 (MES1768 – 3 vials) Lot#: L-NH0010

---

**Cell storage**

- -140°C indefinitely
- Once thawed, useable up to 2 weeks if fed every 4 days

**Product Description**

- Quantity: 3 x 1.0ml (approximately  $7.5 \times 10^6$  cells/cryovial)
- MEF Mitotically inactivated by treatment with mitomycin C
- Both Neomycin and HAT resistant.

**Quality Control**

- All cells were determined to be mycoplasma negative.
- Post-cryostorage viability was determined to be approximately 80%.

**Surface Area Covered**

- $3.5 \times 10^4$  cells/cm<sup>2</sup> to  $9 \times 10^4$  cells/cm<sup>2</sup>